



A fast, easy, and affordable email-based solution to securely send DICOM images and data to physicians through the Internet



Timely and efficient communication is a must in today's healthcare environment. Since 1993 RADinfo SYSTEMS has met increasingly challenging healthcare business needs with affordable products, solutions and services.

PowerPACS DICOMmail allows you to integrate processes across your organization through standard e-mail systems; access and share critical and protected images; facilitate better decision-making and improved clinical outcomes; offer more sophisticated, Internet-enabled, services and gain efficiency to better focus on the delivery of quality patient care.



PowerPACS DICOMmail is RADinfo Systems' DICOM-thru-email solution. DICOMmail comes with our FREE Basic RSVS Viewer and presents the ideal solution for people seeking outside consulting. It allows a user to send DICOM images thru email by simply dragging-and-dropping files; converting JPEG, GIF and/or BMP images into DICOM format, if needed, and then clicking Send. The consulting user can view the images in the Basic RSVS Viewer after a quick and easy Internet download of the DICOMmail application.

PowerPACS DICOMmail was developed to address the need for electronic data exchange within healthcare organizations and their partners, such as associated physicians, consultants, diagnostic imaging centers, clinics, and insurance carriers.

The PowerPACS DICOMmail solution allows the user to anonymize images according to user-defined parameters. The user then sends the images along with a message to healthcare providers, doctors, referring physicians and even patients. When the recipients receive the message they are able to view the image in our free Basic RSVS DICOM Viewer.

PowerPACS[®] DICOMmail



About RADinfo SYSTEMS

Established in 1993, RADinfo SYSTEMS develops and supports Windows-based, DICOM-compliant software for PACS, Teleradiology and image/information systems management. With products installed at hundreds of locations throughout the world, RADinfo SYSTEMS is also a systems integrator, bundling our software with modality and hardware manufacturers' products to provide complete information and image management solutions backed by around-the-clock technical support. In addition to developing the interface software that exists on modalities today, the privately held RADinfo SYSTEMS supplies software products and system development for many major modality and healthcare product vendors.

Features

Digital Distribution

- Supports JPEG, GIF and/or BMP images to send through email
- Flexible solution with Internet application
- Anonymizes images according to user-defined parameters
- An excellent solution to allow referring physicians to view patient images
- Can distribute full-fidelity DICOM images to diagnostic or clinical viewing stations

Extension of the DICOM infrastructure

- Goes beyond medical images
- Maintains structure and quality of each image sent
- Can send video and audio clips with image studies
- Simple drop & drag operation for DICOMmail send
- It has all the advantages that an ordinary e-mail has: light-weighted, flexible and economic
- Handles Structured Reports, ECG and Hemodynamic Waveforms, Encapsulated Documents, Presentation States, Graphic Annotations, Bitmap Plane Overlays, Key Image Notes, etc.



43676 Trade Center Place
Suite 100
Dulles, VA 20166
703-713-3313 phone
703-713-3343 fax
www.radinfosystems.com
info@radinfosystems.com