

PRESS RELEASE FOR IMMEDIATE RELEASE

Media Contact:
Marketing
media@RADinfoSYSTEMS.com
703-713-3313

RADinfo SYSTEMS™ announces the arrival of the PowerPACS® Image Importer, a DICOM importing and management tool for PACS

Dulles, VA, September 16, 2005 – RADinfo SYSTEMS Inc., a leader in affordable RIS/PACS and Teleradiology solutions, announces the launch of its new DICOM importing solution, PowerPACS Image Importer. Healthcare providers now have a set of powerful tools for patient consolidation and management control when importing DICOM medical records from outside facilities into an existing PACS. RADinfo SYSTEMS will demonstrate this technology at booth #2165 during the 2005 Radiological Society of North America (RSNA) conference in Chicago, Nov 27-Dec 2.

PowerPACS Image Importer serves as a virtual supervisor for the importation of diagnostic-quality DICOM images from a CD/DVD as well as non-DICOM media, such as plain film, documents and or other image types, for example, jpegs or tiff files. For seamless consolidation of patient medical records within a PACS, the PowerPACS Image Importer retrieves unique patient order information from a DICOM modality worklist (DMWL) supplied from a HIS/RIS/EMR. Using dynamic matching sequences, these orders are then automatically or manually matched with the newly introduced media. PowerPACS Image Importer automatically updates the DICOM headers; so that the new images will be consistent with records in the existing PACS.

“The PowerPACS Image Importer is an extremely useful set of tools for facilities needing to consolidate patient medical records”, says Dr. Chen-Tai Ma, CEO & President of RADinfo SYSTEMS, “though the DICOM standards committee hasn’t specified a standard way to import DICOM images from a storage media, especially removable media like a CD or DVD, the PowerPACS Image Importer strictly follows the DICOM protocol. The time is near for facilities to have this flexibility to unify their patient records.

About the Image Importer

PowerPACS Image Importer is a low cost, simple solution to import DICOM files from CD/DVDs containing diagnostic quality images, plain film from a film digitizer or jpegs and tiff files from an image directory. The PowerPACS Image Importer connects to any HIS/RIS/EMR, image servers, image workstations and film digitizers to retrieve DICOM patient orders. Either automatically or only on demand, PowerPACS Image Importer uses a smart logic to match the orders with images from the CD or DVD and then import them into the customer’s PACS. The next step for PowerPACS Image Importer is to incorporate a fully scalable PACS migration tools for PACS upgrades from legacy systems. PowerPACS Image Importer runs on any computer with the Microsoft Windows operating system.

About RADinfo SYSTEMS

RADinfo SYSTEMS is leading the revolution towards truly affordable filmless healthcare systems. The company’s products focus on empowering imaging facilities with reliable components of Radiology Information Systems (RIS), Picture Archiving and Communication Systems (PACS), Teleradiology Systems, and Engineering Services at a fraction of the cost of larger competing solutions. RADinfo SYSTEMS develops and integrates DICOM- and HL7-compliant software with other modality and hardware products to deliver affordable solutions to hundreds of small and large healthcare providers worldwide. Professional services include software development, consultation, integration, installation and 24x7 technical support. Delivering Solutions that Meet Your Vision™ since 1993, RADinfo SYSTEMS headquarters are located in Dulles, VA and can be found on the web at www.RADinfoSYSTEMS.com.

PowerPACS® is a federally registered trademark of Radiology Information Systems, Inc. RADinfo Systems™, the RADinfo SYSTEMS circle R logo and “Delivering Solutions that Meet Your Vision™” are trademarks of Radiology Information Systems, Inc.